

## Metallic



Blue



Red



Gold



Violet



Orange



Silver



Green



Steel Blue



Pearl White

## High Gloss



Red



Black



Purple



Yellow



Tan

**Crusader Colors**

These images are a close representation of actual colors, but are not an exact duplication. The colors and fabrics shown are for illustration purposes only, and will vary slightly from the actual finish of the brace.

**Ordering Information:****Part Number****Description**

507100

The Crusader, ACL

507101

The Crusader, CI

507102

The Crusader, PCL

507103

The Crusader, OA

507104

The Crusader, OA &amp; ACL

**Accessories:****Part Number****Description****Size**

50297 X

Under Brace Sleeve

XS - 3X

50622 X

Sports Cover

SM - 2X

590051

Condylar Pad Kit

N/A

590027

0 Degree Flexion Stop

N/A

200757

Knee Brace Bag

N/A

013306

Flexion / Extension Stops

507000

Top Shelf Custom Brace Measuring Instrument



topshelforthopedics.com

1851 East Paradise Road, Suite A • Tracy, CA 95304 • Ph: 866.592.0488 • Fax: 209.834.8832



Crusader™



- Active brace stability with dynamic strapping design neutralizes tibia
- Integrated strapping no external hardware to fail
- Custom fabricated lightweight carbon composite frame provides maximum control & extreme strength
- Specifically designed measuring system captures all needed anatomic measurements allowing superior fit with increased control
- Three lengths available

- Five configurations available - ACL, PCL, Combined Instability, Osteoarthritis, and Combined ACL & Osteoarthritis
- Strap padding has latex free silicone to enhance suspension
- Moisture wicking liner draws moisture away from skin
- Lifetime warranty on brace and hinges
- Fourteen colors available
- Knee Brace Bag

 **TOP SHELF**  
ORTHOPEDICS™

**Our Crusader™ is designed for your crusade!**