

### Crusader Colors

These images are a close representation of actual colors, but are not an exact duplication. The colors and fabrics shown are for illustration purposes only, and will vary slightly from the actual finish of the brace.



### Crusader Ordering Information:

Part Number	Description
507103	The Crusader, OA
507104	The Crusader, OA & ACL

### Knight OA Ordering Information:

Part Number	Description
506302 X	OA,ACL, Short Calf, Left, Medial
506312 X	OA,ACL, Short Calf, Right, Medial
506322 X	OA,ACL, Short Calf, Left, Lateral
506332 X	OA,ACL, Short Calf, Right, Lateral
506342 X	OA,ACL, Standard Calf, Left, Medial
506352 X	OA,ACL, Standard Calf, Right, Medial
506362 X	OA,ACL, Standard Calf, Left, Lateral
506372 X	OA,ACL, Standard Calf, Right, Lateral

Part Number	Description
506382 X	OA, CI, Short Calf, Left, Medial
506392 X	OA, CI, Short Calf, Right, Medial
506402 X	OA, CI, Short Calf, Left, Lateral
506412 X	OA, CI, Short Calf, Right, Lateral
506422 X	OA, CI, Standard Calf, Left, Medial
506432 X	OA, CI, Standard Calf, Right, Medial
506442 X	OA, CI, Standard Calf, Left, Lateral
506452 X	OA, CI, Standard Calf, Right, Lateral

### Accessories:


Part Number	Description	Size
50297 X	Lycra Undersleeve	XS - 3X
590027	0 Degree Flexion Stop	N/A
590052	OA Condylar Pad Kit	N/A
200757	Knee Brace Bag	N/A
013306	Flexion/Extension Stops	N/A
507005	IMI Measuring Instrument	N/A

### The Knight Off the Shelf Sizing Chart

Size #	Size	6" Above Mid Patella	6" Below Mid Patella
X=2	XS	13" - 15"	10" - 12"
X=3	S	15" - 18"	12" - 14"
X=4	M	18" - 21"	14" - 16"
X=5	L	21" - 23"	16" - 18"
X=6	XL	23" - 26"	18" - 20"
X=7	XXL	26" - 29"	20" - 22"
X=8	XXXL	29" - 32"	22" - 24"

Circumference should be taken 6" (15 cm) above mid patella and 6" (15 cm) below mid patella.

# Osteoarthritis Bracing






## Improving your quality of life.....

Osteoarthritis (OA) affects nearly 30 million Americans annually. It is characterized by pain and stiffness of the joints after periods of activity or excessive use. Patients often feel a grating or “catching” during joint movement.

If you have osteoarthritis of the knee, you can take advantage of a wide range of treatment options. The effectiveness of different treatments varies from person to person. The choice of treatment should be a joint decision between you and your physician.

Many people who develop knee pain from OA get good relief from simple lifestyle changes like getting more exercise and from standard treatments like knee braces, shoe inserts, and mild anti-inflammatory medications. Knee braces may be used to limit movement of the knee during recovery from an injury or surgery or relieve knee OA pain by helping correct the alignment and redirecting the body's weight from the affected compartment of the knee joint



### Knight OA Features & Benefits

- Easy to adjust unloading hinge provides an unloading force to the affected compartment
- Integrated strapping system enhances control and reduces profile
- Active brace stability with dynamic strap design provides neutralizing force to the tibia
- Lightweight brace accommodates any time wear
- Graduated thigh reduces medial profile and accommodates bilateral brace wear
- Medial tibial contour controls medial and lateral tibial translation
- Suede moisture wicking liner draws moisture away from skin
- Aerospace-grade aluminum provides rigidity for collision sports
- Two configurations and two lengths are available in ACL & Combined Instability, with short & standard calf
- Latex-free silicone tibial strap pad adds suspension control

### Crusader OA Features & Benefits

- Easy to adjust unloading hinge provides an unloading force to the affected compartment
- Active brace stability with dynamic strapping design neutralizes tibia
- Integrated strapping no external hardware to fail
- Custom fabricated lightweight carbon composite frame provides maximum control & extreme strength
- Specifically designed photo measuring system captures all needed anatomic measurements allowing superior fit with increased control
- Three lengths available
- Available in two configurations ACL / OA, and OA
- Strap padding has latex-free silicone to enhance suspension
- Moisture wicking liner draws moisture away from skin
- Lifetime warranty on brace and hinges
- Fourteen colors available

