

mi-ultra™

Lightweight, dynamic imaging
in the palm of your hand



Innovation that puts physicians first

Trice Medical's continued pursuit to innovate on behalf of physicians takes a leap forward with the NEW **mi-ultra™**. Physicians will have an elegant, lightweight game-changer in their hands, whether diagnosing joint injuries or injecting.

Ease of Use

The ergonomically-friendly mi-ultra™ represents a dynamic solution that is as comfortable to use as it is effective.

Improved Precision

mi-ultra™ is an excellent tool that can provide clarity to physicians whether performing diagnoses or injections.

Versatile Transducers

Available in linear and convex for higher resolution or deeper tissue penetration. Powered by mi-tablet™. Fast image display with no lag.

User-friendly Interface

Easily select desired use and application by using touchscreen display and intuitive icon system.



+ Covered by most major insurance carriers, including Medicare

A Dynamic Imaging Platform

mi-eye® and **mi-ultra™** are a powerful combination of dynamic imaging devices that help physicians diagnose joint injuries with clarity and confidence never seen before.



MI-EYE

A revolutionary handheld arthroscope, mi-eye is a powerful in-office alternative to an MRI.



MI-ULTRA

Compact, lightweight, and easily transportable, mi-ultra is where utility meets high-quality imagery, whether used for diagnosis or injections.



PRODUCT SPECIFICATIONS



mi-ultra™ LR

Linear Ultrasound Transducer

At a mere 80 grams, this linear array transducer rivals many large ultrasound machines. The higher frequency range allows for greater image clarity.

- 1-10 cm Depth Range
- 7.5-15 MHz scan frequency
- Size: 142 mm x 58 mm
- Weight: 80 grams
- 38 mm patient contact



mi-ultra™ CX

Convex Ultrasound Transducer

This lightweight, handheld and portable curved array ultrasound uses a lower frequency permitting deeper tissue penetration.

- 1-21 cm Depth Range
- 3.5-7.5 MHz scan frequency
- Size: 150 mm x 62 mm
- Weight: 110 grams
- 60° Scan Angle
- 49 mm patient contact

With mi-ultra,
treat your patients
with even greater
confidence

Learn more and discover technique guides, videos,
and benefits at tricomedical.com/mi-ultra

Part	Catalog Number
mi-ultra™ LX Linear transducer	T900-020
mi-ultra™ CX Convex transducer	T900-021
mi-tablet	T200-001
95mm mi-eye	T200-095
160mm mi-eye	T200-160
mi-probe, 5-pack	T900-015
mi-eye Procedure Kit	T900-005
<i>Always refer to the products Direction for Use and User Manual for product information and approved use.</i>	