

# Horizon Walkers



**Air Pressure  
Control**

**Adult & Pediatric  
Sizes Available**



**Horizon Walkers are designed to stabilize and protect the lower leg, ankle, and foot. They are easy to apply and comfortable to wear.**

# HORIZON WALKERS

## Order Information

Part #	Description
50170 X	Air Horizon Walker, Tall
50171 X	Air Horizon Walker, Short
50172 X	Horizon Walker, Tall
50173 X	Horizon Walker, Short

Size	Men's Shoe Size	Women's Shoe Size
XS = 2	< 4	< 5.5
SM = 3	< 4.5 - 7	6 - 8
MD = 4	7.5 - 10.5	8.5 - 11.5
LG = 5	10.5 - 12.5	11.5 - 13.5
XL = 6	> 12.5	> 13.5

Part #	Description	Size	Shoe Size
501693	Air Pediatric Horizon Walker	Small	6.5 - 9.5 Child
501694	Air Pediatric Horizon Walker	Medium	Child 10 to Youth I
501695	Air Pediatric Horizon Walker	Large	Youth I - 2.5
501743	Pediatric Horizon Walker	Small	6.5 - 9.5 Child
501744	Pediatric Horizon Walker	Medium	Child 10 to Youth I
501745	Pediatric Horizon Walker	Large	Youth I - 2.5

## Indications

- Stable ankle fracture
- Strains/sprains of the ankle/foot
- Post-surgical stabilization

## Features

- Available in adult sizes with and without air
- High quality molded frame provides for the active patient
- Shock absorbing sole to improve comfort
- Aggressive tread design to prevent slipping when walking
- Low profile rocker bottom sole and wide foot bed to help promote a natural gait
- Integrated air pump for simple inflation and deflation

