

## The Revolution Continues...

The Kodiak fuses the rock-solid dependability patients demand into a compact, lightweight, system that is remarkably quiet and exceptionally easy to use.

With new features, including an ice retainer, lid gasket for reducing spills, and an optional battery pack for portable use, the Kodiak is the result of on-going product innovation.

Through ease-of-use and an enhanced patient experience, the Kodiak leads to improved patient recovery, post surgery.



# The Kodiak... Once again leading the revolution in motorized cold therapy

Splash Proof

Lid gasket

Quiet

Advanced technology for whisper-quiet performance

Simple

Easy-to-use  
Minimal set up required

Comfortable

Ergonomic pads provide exceptional comfort

Reliable

Robust motor for optimal durability

Dependable

Superior filtration system helps prevent clogs

Compact

Lightweight and portable design

Consistent

In-pad flow restrictor maintains consistent temperature



Battery Pack



Ice Retainer  
Retains ice between each ice change



## Optional Battery Pack

With the Kodiak battery pack accessory your patient or athlete can now experience motorized cold therapy in any location, whether it be in the car or on the playing field.

Each pack comes with 4 replaceable AA batteries that will power the Kodiak unit for 10-14 hours depending on application type.



## Intelli-Flo Pads

The POLAR CARE Kodiak features patent exclusive Intelli-Flo pads, providing custom-tailored temperature control specific to each treatment area.

This innovative technology allows for consistent cold delivery, while the ergonomic pad designs provide exceptional comfort and coverage.



For more product or service information, please contact your local sales representative at 1.800.726.9180.

[topshelforthopedics.com](http://topshelforthopedics.com)

POLAR CARE KODIAK	
Part No.	Description
502700	POLAR CARE KODIAK
<b>Combos (Includes 1 Kodiak system and 1 single-patient use pad):</b>	
502701	MULTI-USE COMBO
502702	MULTI-USE COMBO, LONG STEM
502705	SHOULDER COMBO
<b>Individual Pads:</b>	
502710	BACK, INTELLI-FLO PAD
502711	HAND/WRIST, INTELLI-FLO PAD
502712	SHOULDER, INTELLI-FLO PAD
502713	SHOULDER XL, INTELLI-FLO PAD
502714	MULTI-USE, INTELLI-FLO PAD
502715	MULTI-USE, LONG STEM, INTELLI-FLO PAD
502719	ANKLE, COMPRESSION, INTELLI-FLO PAD
502720	3" X 5", INTELLI-FLO PAD
502721	HIP, INTELLI-FLO PAD
<b>Sterile Insulation Barriers (For use only with Intelli-Flo pads):</b>	
502716	KNEE/MULTI-USE/SHOULDER, STERILE POLAR DRESSING
502717	SHOULDER XL, STERILE POLAR DRESSING
502718	ANKLE, STERILE POLAR DRESSING
502722	HIP, STERILE POLAR DRESSING
<b>Accessories</b>	
502890	Y ADAPTER
502891	KODIAK COMPRESSION STRAP
502892	REPLACEMENT TRANSFORMER
502893	KODIAK BATTERY PACK



Distributed by:  
1700 North Chrisman Road  
Tracy, CA 95304  
Phone: 866.323.8395 Fax: 866.323.8612  
[bracing@pacmedical.com](mailto:bracing@pacmedical.com)

1851 East Paradise Road, Suite A  
Tracy, CA 95304  
Ph: 866.592.0488  
Fax: 209.834.8832

