



*Adjustable toe guard
for additional protection.*



Lunar Air Walker

Indications

- Acute ankle sprains
- Achilles tendon repairs
- Post ORIF and/or surgical stabilization
- Stable fractures of the foot & ankle
- Soft tissue injuries of the lower leg
- Stable delayed union or non-union fractures of distal tibia and fibula

Features

- Lightweight nylon composite shell design provides increased stability
- Full shell wraparound for compression
- Independent medial and lateral air bladders control comfort and fit
- Easy access to individual air bladder controls
- Adjustable toe guard for additional protection
- Wide foot bed accommodates swelling and post-op dressing
- Non-skid sole
- Washable liner

Ph: 866.592.0488

Fax: 209.834.8832

1851 E. Paradise Rd., Ste A

Tracy, CA 95304



Lunar Air Walker

Order Information

Part #	Description
501772	Lunar Air Walker, Tall, XS
501773	Lunar Air Walker, Tall, SM
501774	Lunar Air Walker, Tall, MD
501775	Lunar Air Walker, Tall, LG
501776	Lunar Air Walker, Tall, XL

Part #	Description
501782	Lunar Air Walker, Short, XS
501783	Lunar Air Walker, Short, SM
501784	Lunar Air Walker, Short, MD
501785	Lunar Air Walker, Short, LG
501786	Lunar Air Walker, Short, XL

Size	Men's Shoe Size	Women's Shoe Size
XS	Up to 4	Up to 5
SM	4 - 7	5 - 8
MD	7 - 10	8 - 11
LG	10 - 13	11 - 15
XL	13	15

To view our full product line, please visit us at
www.topshelforthopedics.com

