



Experience the Freedom.

Experience **THE REVOLUTION.**



THE REVOLUTION™

Functional Knee Brace

Indications:

- Mild to moderate instability
- ACL, PCL, MCL, LCL injury or surgery
- Hyperextension control

Features:

- Heat formable cuffs - dynamically accommodates leg shape, limits migration
- Lightweight materials; less than 18 oz. - increases patient compliance
- Dynamic strapping system aids in neutralization of the tibia in relation to the femur
- Various flexion and extension levels - accommodates early post-injury or post-operative application
- One year warranty on frame and hinge
- Six months warranty on soft goods padding and straps

The REVOLUTION™ is the newest approach to functional bracing. It is ideally suited to be used in early post-injury or post-operative phase of treatment, following a cruciate or collateral ligament injury.

The REVOLUTION™ is designed to be anatomically contouring, and provides Composite Stability™ allowing the brace and the leg to work as one.



Revolution Knee Brace Ordering

** Measure circumference of thigh 6" (15 cm) above mid patella.*

Left	Size	Right
505402	XS, 13" - 15"	505412
505403	SM, 15" - 18"	505413
505404	MD, 18" - 21"	505414
505405	LG, 21" - 23"	505415
505406	XL, 23" - 26"	505416
505407	2XL, 26" - 29"	505417
Accessories	Size	Description
590019	N/A	Condyle Pad Kit
50297	XS - 2XL	Cotton Lycra Sleeve