



SPD

(Shoulder Positioning Device)



Indications

- Rotator Cuff Repair
- Bankart Repair
- Shoulder Instability
- Total Shoulder Replacement
- Capsule Repair
- Muscle / Ligament Injury
- SLAP Repair
- Shoulder dislocation
- Posterior Capsule Repair
- AC Joint Separation

Features

- Universal lt/rt - reduces inventory
- Malleable waist belt - adjust to all body types
- Internal, neutral & external rotation adjustable to accommodate rehab-specific guidelines
- Adjustable abduction to position the arm during recovery and allows for ventilation to the axilla
- Off-the-neck shoulder strap reduces pain and pressure on the unaffected shoulder
- Adjustable pull tab control for hand grip pad
- Trimable strap-ends allow for adjustment and optimal fit to the wearer
- Washable with mild detergent and air dry
- Allows "Sling Only" configuration to aid in the recovery process when abduction is no longer needed
- Easy one handed buckle closure allows the wearer to don and doff the device without assistance

SPD
(Shoulder Positioning Device)



Order Information

Part Number	Size
502160	Universal

Ph: 866.592.0488


Fax: 209.834.8832

1851 E. Paradise Rd., Ste A

Tracy, CA 95304

To view our full product line, please visit us at
www.topshelforthopedics.com

ARM SLING ASSEMBLY



WAIST ASSEMBLY

