

**Low profile heel height &  
rocker sole for normal gait.**



**Non-Skid Sole**

**Solar Walker**

## Indications

- Soft tissue injuries
- Acute ankle sprains
- Stable fractures of the foot & ankle
- Achilles tendon repair
- Post foot/ankle surgery recovery
- Long term use/chronic ankle/foot pain

## Features

- Malleable aluminum uprights
- XS - XL size form Adult size 4 to size 15
- Removeable & washable liner
- Pneumatic & non-pneumatic liners
- Non-skid rocker bottom sole
- Lower profile height aids in natural gait
- Wide foot bed accommodates swelling and post-op dressing



Air Valve

# Solar Walker



## Order Information

Part #	Description
502682	SOLAR WALKER, TALL, XS
502683	SOLAR WALKER, TALL, SM
502684	SOLAR WALKER, TALL, MD
502685	SOLAR WALKER, TALL, LG
502686	SOLAR WALKER, TALL, XL

Part #	Description
502582	SOLAR AIR WALKER, TALL, XS
502583	SOLAR AIR WALKER, TALL, SM
502584	SOLAR AIR WALKER, TALL, MD
502585	SOLAR AIR WALKER, TALL, LG
502586	SOLAR AIR WALKER, TALL, XL

Part #	Description
502692	SOLAR WALKER, ANKLE, XS
502693	SOLAR WALKER, ANKLE, SM
502694	SOLAR WALKER, ANKLE, MD
502695	SOLAR WALKER, ANKLE, LG
502696	SOLAR WALKER, ANKLE, XL

Part #	Description
502592	SOLAR AIR WALKER, ANKLE, XS
502593	SOLAR AIR WALKER, ANKLE, SM
502594	SOLAR AIR WALKER, ANKLE, MD
502595	SOLAR AIR WALKER, ANKLE, LG
502596	SOLAR AIR WALKER, ANKLE, XL

Ph: 866.592.0488

Fax: 209.834.8832

1851 E. Paradise Rd., Ste A

Tracy, CA 95304

To view our full product line, please visit us at  
[www.topshelforthopedics.com](http://www.topshelforthopedics.com)

Size	Men's Shoe Size	Women's Shoe Size
XS = 2	Up to 4	Up to 5
SM = 3	4 - 7	5 - 8
MD = 4	7 - 10	8 - 11
LG = 5	10 - 13	11 - 15
XL = 6	13+	15+