



Stable Sling 3

Providing

Comfort

Stability

Rotation Control



Indications

- Rotator Cuff Repair
- Bankart Repair
- Shoulder Instability
- Total Shoulder Replacement
- Capsule Repair
- Muscle / Ligament Injury
- SLAP Repair
- Shoulder dislocation
- Posterior Capsule Repair
- AC Joint Separation

Features

- Universal Lt/rt - reduces inventory
- Standard Abduction Pillow, External Rotation Abduction Pillow Options
- Internal, neutral & external rotation adjustable to accommodate rehab specific guidelines
- Abduction to position the arm during recovery and allows for ventilation to the axilla
- Off-the-neck shoulder strap reduces pain and pressure on the unaffected shoulder
- Adjustable hand grip
- Trimmable strap ends allow for adjustment and optimal fit to the wearer
- Washable with mild detergent and air dry
- Allows "Sling Only" configuration to aid in the recovery process when abduction is no longer needed
- Easy one handed buckle closure allows the wearer to don and doff the device without assistance

Stable Sling 3



Part #	Description	Size
502083	Stable Sling 3, Standard Abduction Pillow	Small
502084	Stable Sling 3, Standard Abduction Pillow	Medium
502085	Stable Sling 3, Standard Abduction Pillow	Large
502086	Stable Sling 3, Standard Abduction Pillow	X-Large

Part #	Description	Size
502043	Stable Sling 3, External Rotation Abduction Pillow	Small
502044	Stable Sling 3, External Rotation Abduction Pillow	Medium
502045	Stable Sling 3, External Rotation Abduction Pillow	Large
502046	Stable Sling 3, External Rotation Abduction Pillow	X-Large

Part #	Description	Size
502093	Stable Sling 3 Only	Small
502094	Stable Sling 3 Only	Medium
502095	Stable Sling 3 Only	Large
502096	Stable Sling 3 Only	X-Large

Part #	Description
502060	15 Degree Accessory Wedge

Ph: 866.592.0488
Fax: 209.834.8832
1851 E. Paradise Rd., Ste A
Tracy, CA 95304
www.topshelforthopedics.com



Application Video



External Rotation
Abduction Pillow