



The Protector

Indications:

Post Trans-Tibial amputation protection from injury or impact.

Description:

The ProtectOR is designed to aid in the prevention of injury due to impact or falls following a trans-tibial amputation. It is adjustable in size range to accommodate an individual's anatomy for both the left or right limb. It has been designed to be comfortable and adaptable to a patient's leg with a removable and washable liner. Wearing the ProtectOR will aid in the prevention of flexion contractures and provide a safety barrier over the surgical site when ambulating or transferring from the bed, chair or commode. It allows for easy access and attention to the surgical dressing for visual inspection and changing. No assembly is required once fitted by a practitioner (Circumference can be altered within each size range for 2" of diameter difference).

Order Information

Part #	Description	Size	Circumference
520103	The ProtectOR	Small	15" - 17"
520104	The ProtectOR	Medium	17" - 19"
520105	The ProtectOR	Large	19" - 21"
520106	The ProtectOR	X-Large	21" - 23"

** Measure largest distal circumference of the residual limb.*



Features:

- Easy access for wound care and dressing change
- Adjustable in circumference within size ranges
- Ventilated clamshell design for reduced condensation
- Removable washable liner to reduce contamination and promote hygiene
- Comes with foam inserts to accommodate distal padding
- Semi rigid material allows for flex and contour to limb
- Prevents contracture following surgery
- Optional suspension belt included
- 4 sizes to reduce inventory
- Supracondylar suspension
- No assembly required
- Left – Right Universal
- Latex free



Patent pending

Ph: 866.592.0488

Fax: 209.834.8832

1851 E. Paradise Rd., Ste A

Tracy, CA 95304

To view our full product line, please visit us at

www.topshelforthopedics.com