

Cutter & Bur Guide

Power»Precision»Performance

Call to order:
800.624.4422



Cutters

Actual Size
Hubs Are Color-Coded For Size Only

2.0mm
 2.5mm
 3.0mm
 3.5mm
 4.0mm
 4.5mm
 5.0mm
 5.5mm

Cutter Type	Product No.	Product Description
AGGRESSIVE PLUS	3.5mm 375-534-000 4.0mm 375-544-000 5.0mm 375-554-000 5.5mm 375-564-000	Razor sharp toothed inner blade provides aggressive resection of meniscus, cartilage, and synovium. No teeth on outer window to protect articular cartilage upon entry.
RESECTOR	3.5mm 375-532-000 4.0mm 375-542-000 5.0mm 375-552-000 5.5mm 375-562-000	For general debridement. Particularly effective on articular cartilage. Ideal for removal of periosteum from the acromion. Run in forward for bone removal in acromioplasty and notchplasty.
TOMCAT	3.5mm 375-535-000 4.0mm 375-545-000 5.0mm 375-555-000 5.5mm 375-565-000	The narrow gap between inner and outer cutting surfaces provides aggressive resection of soft tissue in knee and shoulder. Extremely effective on meniscus and plicae.
SCALLOPED	4.0mm 375-546-000 5.5mm 375-566-000	For tough cartilage and ACL stump debridement.
DOUBLE BITE	4.0mm 375-543-000	New, aggressive design shortens time for extensive synovectomy.
SLOTTED WHISKER	4.0mm 375-548-000 5.0mm 375-558-000	For safe, selective debridement of patellar chondromalacia and general selective tissue debridement.
WHISKER CUTTER	4.0mm 375-549-000	Smooths tissue edges.
COUGAR	3.5mm 375-531-000 4.0mm 375-541-000 5.0mm 375-551-000	Open ended design provides both side and end cutting of meniscus, ACL/PCL stump, and plica. Ideal for posterior tears.
END CUTTER	3.5mm 375-737-000 4.0mm 375-747-000 5.0mm 375-757-000	Open ended for meniscal tear debridement.
TRUE END	4.0mm 375-884-000	For direct end cutting in the posterior compartments.

Angled Cutters

Cutter Type	Product No.	Product Description
ANGLED AGGRESSIVE PLUS	4.0mm 380-544-100 4.5mm 380-554-150 5.0mm 380-554-100	Provides rapid tissue removal without clogging. Unique curved design allows increased access to patella, posterior meniscus, and other hard-to-reach areas. No teeth on outer window to protect articular cartilage upon entry.
ANGLED RESECTOR	4.5mm 380-542-150 5.0mm 380-552-100	For general debridement, particularly effective on articular cartilage. Unique curved design allows increased access to patella, posterior meniscus, and other hard-to-reach areas.
ANGLED TOMCAT	4.0mm 380-545-100 4.5mm 380-545-150 5.0mm 380-555-100	Provides rapid, aggressive tissue removal without clogging. Unique curved design allows for general debridement, particularly effective on articular cartilage.
ANGLED DOUBLE BITE	4.5mm 380-543-150	New design provides rapid, aggressive tissue removal without clogging. Unique curved design allows increased access to patella, posterior meniscus, and other hard-to-reach areas. No teeth on outer window to protect articular cartilage upon entry.

Burs

Bur Type	Product No.	Product Description
BARREL BURS (6 flute)	4.0mm 375-941-000 5.0mm 375-951-100 5.5mm 375-951-000	For rapid, aggressive subacromial decompression, notchplasty, and distal clavicle resection without clogging.
BARREL BURS (12 flute)	4.0mm 375-941-012 5.0mm 375-951-112 5.5mm 375-951-012	For smoother bone removal due to increased number of flutes.
ROUND BURS (6 flute)	3.5mm 375-930-000 4.0mm 375-940-000 5.0mm 375-950-100 5.5mm 375-950-000	For rapid, aggressive general abrasion arthroplasty, notchplasty and decompression without clogging.
ROUND BURS (12 flute)	4.0mm 375-940-012 5.0mm 375-950-112 5.5mm 375-950-012	For smoother bone removal due to increased number of flutes.
UNHOODED BARREL BURS (6 flute)	4.0mm 375-941-200 5.5mm 375-951-200	Pre-packaged burs for distal clavicle resection.
UNHOODED ROUND BURS (6 flute)	4.0mm 375-940-200 5.5mm 375-950-200	Pre-packaged burs for distal clavicle resection.
EGG BUR (6 flute)	5.5mm 375-952-000	For rapid, aggressive subacromial decompression and notchplasty.

Speciality Burs

Bur Type	Product No.	Product Description
SLAP™ BUR	4.0mm 375-941-500	For protected use under the labrum while preparing the bony bed of the glenoid in a SLAP lesion.
AUGER BUR	5.0mm 375-450-500	Multi-functional blade resects both bone and soft tissue. Ideal for rapid subacromial decompression and notchplasty when run in forward mode.
SUBCHONDRAL DRILL	4.0mm 375-832-000 NOTE: 1.5mm Drill Tip	Designed for initial treatment of small or unsuspected articular cartilage lesions. Rapidly drills small holes into the subchondral bone promoting vascularity.
LEFT HELIX BARREL BUR (6 flute)	5.5mm 375-953-000	The reverse rotation combined with the left helix flutes draws the bur into the bone during bone resection. This eliminates the bur head from skipping off of the bone and enables surgeons to use the same technique on right and left joints.
LEFT HELIX BARREL BUR (12 flute)	5.5mm 375-953-012	The reverse rotation combined with the left helix flutes draws the bur into the bone during bone resection. This eliminates the bur head from skipping off of the bone and enables surgeons to use the same technique on right and left joints.

Small Joint Cutters & Burs

Cutter/Bur Type	Product No.
AGGRESSIVE CUTTER	2.5mm 375-628-000 3.5mm 375-638-000
FULL RADIUS CUTTER	2.5mm 375-627-000 3.5mm 375-637-000
TOMCAT	3.5mm 375-636-000
HOODED ABRASION BUR	2.0mm 375-641-000 3.0mm 375-647-000

Hip Cutters & Burs

Cutter/Bur Type	Product No.
LONG HIP AGGRESSIVE PLUS	4.0mm 385-544-000
LONG HIP ANGLED AGGRESSIVE PLUS	4.0mm 385-544-100
LONG HIP PEAR BUR (8 flute)	4.0mm 385-943-000