

stryker

Innovation through Integration

- ReconiSense technology allows the CrossFire2 Resection Console to integrate with the CrossFlow Arthroscopy Pump to optimize suction and flow rates for each disposable.
- Connected OR Hub connects seamlessly to EMR, PACS, My Patient Hub iPad App, PowerShare, FTP, and Media Asset Management.
- CrossFire Footswitch has the ability to control SERFAS Energy Probe, Formula Shaver Handpiece, and CrossFlow Arthroscopy Pump features.




Maximized Performance

- CrossBlade Series Shaver Blades reduce particle generation¹, increase resistance to cutter clogging², and maximize cutting performance.³
- CrossFlow's 'Dynamic Mode' allows for more consistent pressure settings to provide optimal distention during resection.
- SERFAS XL Energy Probes are rigid and balanced for hip-joint leverage to ablate and manipulate tissues.



Control Your Environment

- Connected OR Hub allows total control of the entire visualization system from the surgical field through camera head or voice commands.
- Unique surgeon-profiles can be pre-programmed to control features such as: joint pressure & flow, shaver and RF settings, on-screen displays, and footswitch & hand control settings.
- Stryker's ReconiSense TC technology solution is able to monitor the joint temperature state and also provide mitigation measures if necessary.



Portfolio Depth

This portfolio provides a diverse offering of:

- Over 65 different types of shavers and burs
- Suction and non-suction bipolar radio frequency probes
- Inflow only and fully integrated Inflow/Outflow tubing options

	Digital Capture Device
0240200100	CONNECTED OR HUB BASE SYSTEM
	Arthroscopy Consoles
0475100000	CROSSFIRE 2 CONSOLE
0475000100	CROSSFIRE FOOTSWITCH
0375708500	FORMULA 180 SHAVER (HAND CONTROL)
0450000000	CROSSFLOW CONSOLE
	Arthroscopy Disposables
0279351650	SERFAS 50-S SWEEP
0279351653	SERFAS 50-S SWEEP XL
0279401200	SERFAS ENERGY 90-S CRUISE
0475361000	DUAL EDGE SHAVER BLADE - 5.5MM
0475355000	TOMCAT HC SHAVER BLADE - 5.0MM
0475337000	AGGRESSIVE MAX SHAVER BLADE - 3.5MM
0475333000	SMOOTH BITE SHAVER BLADE - 3.5MM
0485860000	XL ROUND DIAMOND BUR 5.5MM
0450000300	CROSSFLOW INTEGRATED CASSETTE TUBING
0279401200	SERFAS ENERGY 90-S CRUISE

Formula series burs

	Standard		Standard
0375941012	STANDARD 12-FLUTE BARREL 4.0MM	0375941200	UNHOODED 6-FLUTE BARREL 4.0MM
0375951112	STANDARD 12-FLUTE BARREL 5.0MM	0375951200	UNHOODED 6-FLUTE BARREL 5.50MM
0375951012	STANDARD 12-FLUTE BARREL 5.5MM	0375940200	UNHOODED 6-FLUTE BARREL 4.0MM
0375940012	STANDARD 12-FLUTE BARREL 4.0MM	0375950200	UNHOODED 6-FLUTE BARREL 5.5MM
0375950112	STANDARD 12-FLUTE BARREL 5.0MM	0375952000	EGG 6-FLUTE 5.5MM
0375950012	STANDARD 12-FLUTE BARREL 5.5MM		
0375941000	AGGRESSIVE 6-FLUTE BARREL 4.0MM		
0375951100	AGGRESSIVE 6-FLUTE BARREL 5.0MM		
0375951000	AGGRESSIVE 6-FLUTE BARREL 5.5MM		
0375930000	AGGRESSIVE 6-FLUTE BARREL 3.5MM		
0375940000	AGGRESSIVE 6-FLUTE BARREL 4.0MM		
0375950100	AGGRESSIVE 6-FLUTE BARREL 5.0MM		
0375950000	AGGRESSIVE 6-FLUTE BARREL 5.5MM		

	Specialty
0375941500	SLAP 6-FLUTE 4.0MM
0375450500	AUGER 5.0MM
0375832000	SUBCHONDRAL DRILL 4.0MM
0375953000	LEFT HELIX 6-FLUTE BARREL BUR 5.5MM
0375953012	LEFT HELIX 12-FLUTE BARREL BUR 5.5MM

TPS small-joint series cutters and burs

	Cutters and burs
0275627000	2.5MM SMALL-JOINT FULL RADIUS CUTTER
0275637000	3.5MM SMALL-JOINT FULL RADIUS CUTTER
0275628000	2.5MM SMALL-JOINT AGGRESSIVE CUTTER
0275638000	3.5MM SMALL-JOINT AGGRESSIVE CUTTER
0275641000	2.0MM SMALL-JOINT HOODED ABRASION BUR
0275647000	3.0MM SMALL-JOINT HOODED ABRASION BUR
	Cutter and bur hardware
0242100000	SMALL-JOINT INSTRUMENTATION KIT
0275035302	2 3.5MM SMALL-JOINT CUTTER CANNULA LOCKS IN 2.3MM AND 2.7MM SMALL-JOINT SCOPE AND ALL 275 SERIES SMALL-JOINT CUTTERS AND BURS (EXCEPT 3.0MM HOODED ABRASION BUR)
0275035100	3.5MM TROCER
0275035200	3.5MM OBTURATOR
0275035300	SEALING CAP
0279250101	2.5MM MICRO-CLAW SERFAS ENERGY SMALL-JOINT PROBE DESIGNED FOR DIRECT END CUTTING AND ABLATION IN SMALL JOINTS
0279250201	2.5MM MICRO-BRUSH SERFAS ENERGY SMALL-JOINT PROBE
0279251101	AN EFFECTIVE ALL-PURPOSE SMALL-JOINT PROBE 2.5MM AARDVARK SERFAS ENERGY SUCTION PROBE BENDABLE SUCTION PROBE FOR TIGHT JOINT SPACES

Cross blade series cutters

	Premium
0475331000	CROSSBLADE DUAL EDGE 3.5MM
0475341000	CROSSBLADE DUAL EDGE 4.0MM
0475355000	CROSSBLADE TOMCAT HC 4.0MM
0475355000	CROSSBLADE TOMCAT HC 5.0MM
0475361000	CROSSBLADE DUAL EDGE 5.5MM
0475345000	CROSSBLADE AGGRESSIVE MAX
0475333000	CROSSBLADE SMOOTH BITE
0475541000	CROSSBLADE BURSECTOR 4.0MM
0475561000	CROSSBLADE BURSECTOR 6.0MM
0485840000	CROSSBLADE DIAMOND BUR XL 4.0MM
0485860000	CROSSBLADE DIAMOND BUR XL 5.5MM

Formula series XL cutters and burs

	Hip length
0385544000	AGGRESSIVE PLUS, LONG HIP CUTTER 4.0MM
0385544100	ANGLED AGGRESSIVE PLUS, LONG HIP CUTTER 4.0MM
0385545000	TOMCAT, LONG HIP CUTTER 4.0MM
0385545100	ANGLED TOMCAT, LONG HIP CUTTER 4.0MM
0385552000	STRAIGHT RESECTOR, LONG HIP CUTTER 4.0MM
0385943000	LONG HIP 8-FLUTE PEAR BUR 4.0MM
0385953008	LONG HIP 8-FLUTE PEAR BUR 5.5MM
0385950008	LONG HIP 8-FLUTE ROUND BUR 5.5MM
0385951008	LONG HIP 8-FLUTE BARREL BUR 5.5MM

Cutter cannula set

	Standard
0270713100	TROCER SE-4
0270713200	OBTURATOR SE-4
0270713300	CUTTER CANNULA WITH HOLE 5MM
0272099000	CANNULA ADAPTER COMPETITIVE
0272100100	DYONICS ADAPTOR BRIDGE, ONE ROTATING STOPCOCK
0272100200	DYONICS ADAPTOR BRIDGE, TWO ROTATING STOPCOCKS
0275035100	SMALL JOINT TROCER 3.5MM
0275035200	SMALL JOINT OBTURATOR 3.5MM
0275035300	SMALL JOINT CAP
0275035302	SMALL JOINT CANNULA
0275719000	CUTTER/BUR CANNULA 5.6MM
0275719100	TROCER 5.6MM
0275719200	OBTURATOR 5.6MM
0275719300	CUTTER/BUR CANNULA (WITH HOLE) 5.6MM
0275720000	BUR CANNULA 5.0MM
0275720100	TROCER 5.0MM
0275720200	OBTURATOR 5.0MM
0275721000	BUR CANNULA 6.0MM
0275721100	TROCER 6.0MM
0275721200	OBTURATOR 6.0MM
0275722000	CUTTER CANNULA 4.5MM
0275722100	TROCER 4.5MM
0275722200	OBTURATOR 4.5MM
0275722300	CUTTER CANNULA WITH HOLE 4.5MM
0333000001	PACKAGE SUCTION/IRRIGATION ADAPTER WITH STOPCOCK
0333070000	70MM IN/OUTFLOW CANNULA KIT
0333070001	DRAIN CANNULA 70MM
0333070002	TROCER FOR 70MM DRAIN CANNULA
0333070003	OBTURATOR FOR 70MM DRAIN CANNULA
0333100000	10MM IN/OUTFLOW CANNULA KIT
0333100001	10MM DRAIN CANNULA
0333100002	TROCER FOR 10MM DRAIN CANNULA
0333100003	OBTURATOR FOR 10MM DRAIN CANNULA
0342010133	FLAT WASHER
0375035100	3.5MM FORMULA SERIES TROCER
0375035302	3.5MM FORMULA SERIES, SMALL JOINT CUTTER CANNULA
0272100000	DYONICS FLUID SEAL ADAPTOR

Formula series cutters

	Standard
0375532000	RESECTOR 3.5MM
0375542000	RESECTOR 4.0MM
0375552000	RESECTOR 5.0MM
0375562000	RESECTOR 5.5MM
0375534000	AGGRESSIVE PLUS 3.5MM
0375544000	AGGRESSIVE PLUS 4.0MM
0375554000	AGGRESSIVE PLUS 5.0MM
0375564000	AGGRESSIVE PLUS 5.5MM
0375535000	TOMCAT 3.5MM
0375545000	TOMCAT 4.0MM
0375555000	TOMCAT 5.0MM
0375565000	TOMCAT 5.5MM
0375546000	SCALLOPED 4.0MM
0375541000	COUGAR END 4.0MM
0375747000	END 4.0MM
0375543000	DOUBLE BITE 4.0MM
0375548000	SLOTTED WHISKER 4.0MM
0375549000	WHISKER 4.0MM

	Angled
0380542150	RESECTOR 4.5MM
0380544100	RESECTOR PLUS 4.0MM
0380544150	RESECTOR PLUS 4.5.0MM
0380554100	RESECTOR PLUS 5.0MM
0380545100	TOMCAT 4.0MM
0380545150	TOMCAT 4.5MM
0380555100	TOMCAT 5.0MM
0380543150	DOUBLE BITE 4.5MM

Literature number:1000903377 Rev.A

1. Fulcrum Particle Generation Test: TP11189 Data: DHD12753; TR14605; TR16038; DHD1310
2. Dual Mode Clogging Test: TP10415 Data: DHD13100; TR16038
3. Ligament Consumption Test: TP11560 Data: TR15927; DHFD11511; DHD13100; TR16038